



Instruction for Heartscope & Sleepscope electronic referral templates for GP's



Best Practice
An evolution in medical software

Instructions

- Step 1** Open Patient's file & create referral as usual, Select "Blank Document"
- Step 2** Once blank document opens, select the "Blank Document" with writing
- Step 3** The "Word Processor Templates" will open - Scroll down and search for which Referral you would like: (**Heartscope Cardiology Referral** or **Sleepscope Sleep and Respiratory Referral**)

Healthlink EDI: heartvic

To make a booking please call: **1800 202 111**

fax referral to 8669 4575

For more information about our services, to request referral set up for Best Practice or Medical Director or for printed referral pads, please contact **Amanda Spark** on **7002 1480**